



MULTI ELEMENT PROBES

The high throughput multi-element probes increase productivity and minimize repetitive tasks by processing numerous samples simultaneously. Units are available with 4, 8, 16 and 24 tips and are compatible with either the 500 or 750 watt systems. Custom formatted multi-element probes are available upon request. The replacement tips are Part no. 630-0660.

4-ELEMENT



PART NO.
630-0559

8-ELEMENT



PART NO.
630-0586

16-ELEMENT



PART NO.
630-0699

24-ELEMENT



PART NO.
630-0579

TIP DIAMETER
1/8" (3mm)

VOLUME
0.5-15mL

AMPLITUDE
120µm

TIP SPACE
18mm



STAND PART NO.
830-00130

HEAVY DUTY SUPPORT ASSEMBLY

Supports the converter and multi-element probe with minimum deflection. Includes lab jack. Part no. 830-00130 is recommended when working with any multi-element horn.